


## CmpSpots - 2D Spot Pairing



<http://open2dprot.sourceforge.net/CmpSpots>

Revised: 12-04-2004, P. Lemkin, NCI-Frederick

# Outline

1. Overview
2. Algorithm
  - 2.1 Landmark sets
  - 2.2 Pairing features dP and dL
  - 2.3 Primary spot pairing
  - 2.4 Secondary pairing
3. Web page
4. Report window GUI
5. Options Wizard GUI
6. Image Viewer - GUI
7. Image Viewer - Examples

## Overview

- CmpSpots is an open-source Java 2D quantified spot pairing program for matching spots between two lists of spots from 2D data samples.
- “Spot” data could be from 2D electrophoretic gel images, 2D LC-MS peptide clusters, or compatible data source, etc.
- There may be different numbers of spots in the two samples implying missing data or N:1 or 1:N spot matching.
- CmpSpots is a step [4] module in the pipeline analysis for the Open2Dprot project.
- Spot data may be used from the Open2Dprot Seg2Dgel 2D PAGE gel spot quantification program or compatible source.

## Overview - Continued (2)

- Input spot lists are read from Sample Spot-list Files (SSF) XML (or tab-delimited) files and paired-spot lists are generated in XML, tab-delimited, or human readable formats as Sample Comparison Files (SCF).
- Uses a sample description experiment “accession” database listing the gels to be used and their region of interests. This DB will be replaced with a MIAPE compliant XML database format and a tool added to Open2Dprot.
- Uses a set of predefined landmarks between the sample gel and a reference gel from a landmark database. This DB will be interactively defined as an XML DB in the future using a new Open2Dprot tool.

### Overview - Continued (3)

- The program may be run either interactively (-gui) with a graphical user interface (GUI) or under an OS shell command to implement batch (-nogui).
- Run-time options are specified as Unix-style command line '-' prefixed input switches so the program is also usable under batch. These may be set by an Options Wizard window.
- In the GUI version, after the spot pairing is finished, the user has the option of interactively viewing the paired spot data generated by the spot pairing with the Image Viewer window.
- Web site contains documentation, downloads and examples.

5


### Algorithm (1)

1. Three parameters are required: thrSP and thrPP the maximum distance between warped putative spots for them to be labelled as Sure-Pairs and as Possible-Pairs; number of alternative landmarks to check when doing secondary pairing.
2. Two input spot lists are read from XML files for a Reference sample and the sample being matched.
3. The landmarks for these samples are read from a landmark database and landmarks effective radii computed.
$$R_k = 1/2 \text{ Min}(\text{distance}(\text{Lm}_k, \text{Lm}_m))$$
for all landmarks k,m where k not m
4. Spots are assigned to closest landmark as landmark sets, as well as noting their next-nearest landmarks.

6

## Effective Radii of Landmarks

Estimates of spot pairs within landmark radii are more reliable




7

## Definition of dP - distance between spots in a pair

The distance between paired spots is  $dP_{kij}$  and is computed after the two spots  $(s_{1i}, s_{2j})$  are mapped to  $(s'_{1i}, s'_{2j})$  in the same sample space relative to landmark k.

$$dP_{kij} = |s'_{1i} - s'_{2j}|.$$



8


### Definition of dL - distance from a pair to it's landmark

The distance between a spot-pair and the landmark is  $dL_{kij}$  and is computed after the two spots ( $s_{1i}, s_{2j}$ ) are mapped to ( $s'_{1i}, s'_{2j}$ ) in the same sample space relative to landmark k.

$$dL_{k1i} = \text{distance}(\text{LMk}, s'_{1i}),$$

$$dL_{k2j} = \text{distance}(\text{LMk}, s'_{2j}),$$

$$dL_{kij} = \min(dL_{k1i}, dL_{k2j}).$$



9

### Algorithm (continued) (2)

5. For each landmark, spots are optimally paired between the corresponding landmark sets:

a) mutually-paired spots  $i, j$  assigned sure-pair (SP) labels if  $(dP_{kij} \leq \text{thrSP})$  and  $(dL_{kij} \leq R_k)$ .

b) mutually-paired spots assigned possible-pair (PP) labels if  $(dP_{kij} \leq \text{thrPP})$  or  $((dP_{kij} \leq \text{thrSP}) \text{ and } (dL_{kij} > R_k))$ .

c) N:1 or 1:N paired spots assigned ambiguous-pair (AP) labels if N spots in one sample pair with 1 spot in other sample.

d) Unresolved Spots (US) are what remains unpaired.

10


## Primary Pairing Cases

### Case

- [1] US is unresolved spot (no dP),
- [2] SP is  $dL \leq R_k$  and  $dP \leq dTsp$ ,
- [3] PP is  $dL \geq R_k$  and  $dP \leq dTpp$ ,
- [4] PP is  $dL \leq R_k$  and  $dP \geq dTsp$  and  $dP \leq dTpp$ ,
- [5] PP is  $dL \geq R_k$  and  $dP \leq dTsp$ ,
- [6] PP is  $dL \geq R_k$  and  $dP \leq dTpp$ . For the other spot AP' is  $dL' \geq R_k$  and  $dP' \leq dTpp$  and  $dP' \geq dP$

11

## Primary Pairing Cases



12

### Algorithm (continued) (3)

6. Secondary spot pairing can be used to further resolve AP and US labels in adjacent landmark sets into SP or PP labels which are then placed in either of the two sets. There are four cases: (a) two unresolved spots (US and US), (b) two ambiguous pairs (AP and AP), (c-d) one ambiguous spot (AP) and one unresolved spot (US).

If the new spot pair is better than its previous pairing, it is upgraded to a SP or PP pairing label and is put into whichever landmark set has the smallest dL for the putative pair.

7. The data is then saved in a SPF output file for further processing by other programs such as the Open2Dprot database builder.


13

### Secondary Pairing Cases

- A second pass is optimized the spot pairing to further resolve AP and US labels in adjacent landmark sets into PP labels that are then placed in either of the two sets.
- There are four cases (see Reference Manual for details):
  - (a) two unresolved spots (US and US),
  - (b) two ambiguous pairs (AP and AP),
  - (c-d) one ambiguous spot (AP) and one unresolved spot (US).
- The new spot pair is upgraded to a SP or PP pairing label and is put into whichever landmark set has the smallest dL for the putative pair.
- Each landmark set is tested. Each spot has a list of adjacent landmarks, so these are tested to see if pairing can be optimized for that spot.

14

## Secondary Pairing Cases



15

<http://open2dprot.sourceforge.net/CmpSpots>

**Table of Contents**

**CmpSpots**

- Home
- Introduction
- Overview (PDF)
- Contributors

**Reference manual**

- Licenses
- PDFs
- Tutorials
- Demos & screen shots
- Version
- Downloads
- Status
- Revision history
- FAQ
- References

**Source code**

- Project summary
- CVS access
- Browse CVS
- Javadocs
- Files mirror

**Other Web resources**

Hosted at Open2Dprot  
[open2dprot.sourceforge.net](http://open2dprot.sourceforge.net)

Powered by

[Contact us](#)

Revised: 08/05/2004

**CmpSpot - 2D Database Spot Pairing**

**CmpSpots - 2D Database Spot Pairing**  
pre-alpha version

Welcome To CmpSpots

<http://www.lecb.ncifcrf.gov/Open2Dprot/CmpSpots/>  
and  
[FUTURE] <http://open2dprot.sourceforge.net/CmpSpots>

**In Table of Contents:**

- \* [Introduction](#)
- \* [Overview \(PDF\)](#)
- \* [Reference Manual](#)
- \* [Demos & Screen shots](#)
- \* [Download](#)
- \* [Source code](#)

16




## Downloading and Installing CmpSpots

- You may download the pre-alpha version and install it on your computer from the Web site.
- Currently, CmpSpots is hardwired to start with the demo gel and with the -gui switch. However, you can override this with the **Edit options** popup wizard.

17

## Installing CmpSpots on your computer



18

## CmpSpots Report Window Interface in -gui mode

The screenshot shows the CmpSpots V0.0.4 pre-Alpha GUI window. The title bar indicates the date 2004/09/05 10:53:59 and the revision 1.3. The menu bar includes File, Edit, View, and Help. The main text area contains a detailed report of image file processing, including file names, sizes, and pairing statistics. The bottom of the window features a set of command buttons: Done, Clear, SaveAs, Stop spot pairing, Edit options, Pair spots, Image viewer, and Close. Annotations with arrows point to the 'Pull-down menus' (File, Edit, View, Help), the 'Report window text' (the main report area), and the 'Command buttons' (the bottom bar).

**Pull-down menus**

**Report window text**

**Command buttons**

19

## CmpSpots Report Window Interface in -gui mode

This screenshot is identical to the one above, showing the CmpSpots V0.0.4 pre-Alpha GUI window. It displays the same report of image file processing and pairing statistics. The annotations highlight the 'Pull-down menus' (File, Edit, View, Help), the 'Report window text' (the main report area), and the 'Command buttons' (the bottom bar: Done, Clear, SaveAs, Stop spot pairing, Edit options, Pair spots, Image viewer, Close). The page number 20 is visible in the bottom right corner.

**Pull-down menus**

**Report window text**


**Command buttons**

20

## Option Wizard to Edit Command Options

- The user may modify the command line switch options using the Option Wizard by pressing the **Edit Options** button.
- You may save the new options in a "*CmpSpots.properties*" file in the current project directory .
- The new options are then used as the new default switch options in subsequent running of CmpSpots.

21



22

## Image Viewer Window

- Interactively view the spot pairing results overlayed on the original images using the Image Viewer.
- Select spot to examine the features of the spot-pair.
- View (landmarks, pairing-labels, spot-pairing vectors, landmark radii) as overlays.
- View subsets of above for currently selected landmark.
- Dynamic histograms spot-pairing statistics (dP, dL, pairing-labels, LM set size) features.
- Filter by pairing-labels and/or histogram features.

23

## Data used in the following examples

**Rsample** (G1=Leukemia-AML\*) has 933 in all landmark sets,  
**Sample** (G2=Leukemia-ALL\*) 2143 spots in all landmark sets.

After Initial pairing:

US 853  
SP 242  
PP 1028  
AP 792  
CP 0  
EP 0  
 $0.5(SP+PP)/(|G1| \text{ MIN } |G2|)=68.06\%$

After secondary pairing:


US 834  
SP 242  
PP 1082  
AP 757  
CP 0  
EP 0

After secondary pairing:  $0.5(SP+PP)/(|G1| \text{ MIN } |G2|)=71.0\%$   
mean dP(SP+PP)=4.66, mean dP'((|G1|+|G2|)/(SP+PP))=8.39

\*Data from CmpSpots demo database


24

## Image Viewer - no overlays, no spot filters




25

## Image Viewer - landmarks




26


## Image Viewer - landmarks and their effective radii




## Image Viewer - labels overlay, no spot filters



## Image Viewer - labels and landmarks radii overlays




## Image Viewer - labels overlay, current landmark set "P"






## Image Viewer - labels & radius overlay, landmark set "P"



31


## Image Viewer - vector difference overlay, no filter




32




## Image Viewer - vector overlay, landmark set "V"



## Image Viewer - labels & radius overlay, landmark set "V"




## Image Viewer - landmark sets "A, B, C, D"




35

## Image Viewer - histogram of pairing features




36

## Image Viewer - spot-pairs filtered by dP, dist. Between spots



## Image Viewer - spot-pairs filtered by dL, dist. LMs to pairs



## Summary

- CmpSpots is an open-source 2D spot pairing Java program freely available at <http://open2dprot.sourceforge.net/CmpSpots>
- Useful for pairing “spots” in 2D samples (such as 2D PAGE gels, 2D LC-MS and other images with similar types of data).
- It may be used as one of the step [4] alternative modules in the analysis pipeline in the Open2Dprot project at <http://open2dprot.sourceforge.net>